

Volatile Chemicals and Class II Type A2 “Recirculated” BSCs: How Much is Safe?

The logo consists of the word "BAKER" in a bold, white, sans-serif font, centered within a blue rounded rectangle.

BAKER

Environments For Science™

Kara Held, Ph.D., Science Director

Kara F. Held, Ph.D.

- **B.S.** in Biomedical Sciences
- **Ph.D.** in Cell and Molecular Biology, focus in Pharmacology
- **Postdoctoral** Training in Vascular Biology and Therapeutics
- **Lab Manager/Researcher/Lab Safety Officer** in Stem Cell and Regenerative Biology
- kheld@bakerco.com

QUINNIPIAC
UNIVERSITY



The
UNIVERSITY
of VERMONT

Yale University



HARVARD
UNIVERSITY

HSCI

HARVARD STEM CELL
INSTITUTE®

Presentation agenda

- Types of BSCs and HEPA filters
- Volatile Chemical use in Biology
- Calculations of Airflow and Chemical Concentrations
- How much can I use?
- Recirculated A2 vs. 100% exhausted B2
- NSF safety standards – “minute” amounts


Biosafety Cabinets (BSCs)

- A ventilated enclosure for work with biohazard agents assigned to biosafety levels 1 through 4, as per the BMBL.
- Provides levels of **CONTAINMENT**
 - Personnel protection
 - Product protection
 - Environmental protection




HEPA filtration

- **High Efficiency Particulate Air**
 - Filter out 99.97% (or better) of 0.3 μm mass median particles of DOP (Dispersed Oil Particulate)




IEST-RPC **CANNOT filter out gasses and vapors!**

Biosafety Cabinet Classification



Class II Type A2 BSCs


- Personnel, Product and, Environmental Protection
- Minimum 100 fpm intake air
- Partial recirculation



Volatile Chemicals in an A2

- If a volatile chemical is introduced in the work area, all air becomes contaminated.
- How much is safe?
 - Vent to room vs. Canopy
 - Explosion limits for BSC
 - Chemical exposure limits
 - “Minute” amounts per NIOS Standard 40

Canopy/“Thimble”



Why can't we just assign a value for a maximum amount of chemical allowed?

- NSF wanted to put a 250 mL limit on volatile chemicals for A2s.
- Not all chemicals are created equal!
 - Different harmful levels
 - Different evaporation rates
- Depends on the surface area of the container used
 - Spill vs. using a capped bottle?
- What about new chemicals not yet evaluated?
- Depends on the amount of air being diluted into


Calculating Airflow

- Airflow is designated by “Q”
- Airflow is the product of velocity and area
 - $Q = V * A$
 - $V = 100$ fpm intake velocity,
 - $A =$ front access opening
(e.g. 8”sash, 4’ cabinet = 2.67 ft²)
 - $Q = 267$ cfm



*Round numbers are used for simplicity. Actual values may vary depending on manufacturer.


Mathematically Derived and Empirically Proven Equations



- Class II Type A cabinets
 - $C_d = C_i + \frac{X}{Q_i}$
 - If C_i in lab = 0, then $C_d = \frac{X}{Q_i}$
- Example:
 - $C_d = \frac{1558 \text{ mg/min}}{7.1 \text{ m}^3/\text{min}}$
 - $C_d = 219 \text{ mg/m}^3$
- Confirmed with toluene empirically

Calculating airflows and volumes

- The amount of chemical in the air is determined by mass flow rate (\dot{m})
 - $\dot{m} = Q * C$
 - Mass flow rate = airflow*concentration
 - Example: $\dot{m}_i = Q_i C_i$



Mass Flow Rates


- Solved as an infinite series of recirculations:

- $\dot{m}_t = \dot{m}_d + \dot{m}_i + \dot{m}_x$

- Where $\dot{m}_{d_n} = R\dot{m}_{t_n}$

- Thus:

- $\dot{m}_{d_\infty} = \left(\frac{R}{1-R}\right) (\dot{m}_i + \dot{m}_x)$



Mass Flow Rate Solution

- Using the empirical data from the toluene example:
- $\dot{m}_{d_{\infty}} = \left(\frac{R}{1-R} \right) (\dot{m}_i + \dot{m}_x)$
- Toluene was released at 1.8 mL/min = 1558 mg/min; 70% recirculation
- $Q_d C_d = 2.33(Q_i C_i + \dot{m}_x)$
- $16.0 \text{ m}^3/\text{min}(C_d) = 2.33(0 + 1558 \text{ mg/min})$
- $C_d = 226 \text{ mg/m}^3$


219 mg/m³ (empirically confirmed) \approx **226 mg/m³** (calculated)

How MUCH can I use?

- Dilution Ventilation equation:

- $$ER = \frac{Q_i * MW * LEL * 473}{403 * SG * S_F * 100}$$

- S_F generally = 4 for 25% of LEL,
more cautious users $S_F = 10$
(10% of LEL)



How much can I use?

- Dilution Ventilation equation:

- $$ER = \frac{Q_i * MW * LEL * 473}{403 * SG * S_F * 100}$$

- *Example: isopropanol in 4ft A2 BSC*

- $$ER = \frac{267 \text{ cfm} * 60 \text{ g/mol} * 2.0\% * 473}{403 * 0.786 * 4 * 100}$$

- **$ER = 120 \text{ mL/min}$**

Calculating a spill

100 mL of 70% isopropanol spilled on the work surface of an 6' A2 BSC.

- Covered 0.313 m² area
- Took 1 hour to evaporate



250 mL of 70% isopropanol covering the entire the work surface of a 6' A2.

- Covers 0.79 m² area
- Will take 1 hour to evaporate


Am I safe?

As Calculated

120 mL/min 100%
isopropanol

-or-

171 mL/min 70% isopropanol



As Spilled

250 mL/hr 70% isopropanol

-or-

4.2 mL/min 70% isopropanol

**40-fold less than
the safe limits!**

How long does it take to clear?

- Reduce contaminated air to a nominal or background amount.
- One recirculation takes ~4 seconds*
- Helium injected to 73,000 ppm
- Measured to background = **60 to 90 seconds**
- $C_X * R^n = C_{X,low}$
- $73,000 \text{ ppm} * (0.7)^n = 10 \text{ ppm}$
- $n = 25 \text{ recirculations}$
- $25 * 4s = \mathbf{100 \text{ seconds}}$

*Round numbers are used for simplicity. Actual values may vary depending on manufacturer.
Empirical study done by The Baker Company

Is a 100% exhausted Type B2 cabinet really needed for volatile chemical use?

- B2 has greater chance of **failure**, because it is reliant on the building's exhaust system.
 - A2 can vent to room or through canopy connection
- B2 has a greater **operating** and **installing cost** relative to A2
- B2 is more **complex** to operate properly
 - Room pressurization, room size, alarms, filter loading lessening static pressure
- Do a risk assessment to determine which is right for your needs!

Determination of “Minute”

- The least amount of chemical needed for your experiment currently being conducted.
- Below explosion and exposure limits
- In appropriate containers, preferably capped.
- ALWAYS perform a Risk Assessment

Acknowledgments

- Gary Hazard, Baker Director of Manufacturing
- David Eagleson, Baker President
- Dan Ghidoni, NESA Sales Manager
- Dan Eagleson, Baker Senior Vice President
- Aaron Johnson, Baker Product Development Supervisor
- David Stuart, Ph.D., Eagleson Institute
- Allison Weaver, Baker Student Intern

References

- Stuart, D.G., First, M.W., Jones, Jr., R.L., and Eagleson, J.M. (1983) Comparison of Chemical Vapor Handling by Three Types of Class II Biological Safety Cabinets. Particulate & Microbial Control
- Caravanos, J. (1991) Quantitative Industrial Hygiene: A Formula Workbook. Am Conf of Gov Ind Hygienists

BAKER

Environments For Science™

Questions?