



UNITED STATES ARMY MEDICAL RESEARCH INSTITUTE OF INFECTIOUS DISEASES

Biodefense solutions to protect our nation

Summary of Validated and Verified Viral Inactivation Methods

Dr. David Harbourt
October 16, 2017



USAMRIID

United States Army
Medical Research Institute
of Infectious Diseases

Biodefense solutions to protect our nation

David Harbourt, PhD, CBSP, SM(NRCM), RBP
Biosafety Officer
David.e.harbourt.civ@mail.mil
301-619-9875
www.usamriid.army.mil



Disclaimer

Biodefense solutions to protect our nation

The views, opinions and findings contained herein are those of the author and should not be construed as an official Department of the Army position, policy or decision unless so designated by other documentation.



Outline

Biodefense solutions to protect our nation

- Background on requirements
- Buffer AVL Validation
- SDS Buffer Validation
- Formalin Validation
- Trizol LS Validation
- Conclusions



Background

Biodefense solutions to protect our nation

- In July 2016, Army Directive 2016-24 was published
 - Required validation of any BSAT inactivation process with in house testing data
 - Any inactivated material movement from registered to unregistered space could not occur without Biosafety, RO and Command approval
- Many requirements in AD 2016-24 are now found in latest 42 CFR updated regulations
- USAMRIID needed to develop validation testing data for inactivation protocols to continue operations



Background

Biodefense solutions to protect our nation

- Many impediments to developing validation datasets
 - Time
 - Resources (monetary and personnel)
 - Methods to remove toxic chemicals affecting cell culture
 - Selection of viruses to serve as surrogates



Background

Biodefense solutions to protect our nation

- Three methods to remove toxic chemicals for virus cell culture
 - Concentrating columns
 - Desalting columns
 - Dilution

VALIDATION STUDIES



Buffer AVL Methodology



Biodefense solutions to protect our nation

- TC83, a RG2 enveloped virus, was grown and concentrated with an Amicon Ultra-15 Centrifugal Filter Unit with a molecular weight cutoff size of 100kDa
- TC83 selected for high viral titer potential and quick propagation relative to other viruses
- Non-concentrated viral titer was 9.8E9 pfu/mL
- Concentrated viral titer was 4.4E11 pfu/mL
- NHP plasma was spiked with the concentrated virus at a 3:1 ratio to create plasma containing 1.1E11 pfu/mL



Buffer AVL Methodology

Biodefense solutions to protect our nation





Buffer AVL Methodology

Biodefense solutions to protect our nation


- Supernatant was collected from the dishes and applied to Hela cells in a 96 well plate format for 18 hours before fixing in 10% neutral buffered formalin for one hour.
- Formalin was washed from the plates twice with PBS
- Plates were blocked with 3% BSA in PBS for 1 hour at ambient temperature
- Plates were incubated with mouse 1A4A anti-VEE antibody at a 1:1000 dilution in 3% BSA in PBS for 1 hour at ambient temperature
- Plates were washed 3x with PBS
- Secondary antibody of DyLight goat anti-mouse 488 (ThermoFisher catalog #35502) was added to the plates at a 1:1000 dilution in 3% BSA in PBS and incubated for 1 hour at ambient temperature
- Plates were washed 3x with PBS and Hoescht and Cell Mask Deep Red were added at a 1:10000 dilution in PBS to stain cell nuclei and cytoplasm
- Image analysis was obtained using an Opera confocal reader (model 5025-Quadruple Excitation High Sensitivity [QEHS])



Buffer AVL Results

Biodefense solutions to protect our nation

Percent Positive Cells after Passaging Inactivated Samples






Buffer AVL Results

Biodefense solutions to protect our nation

Cell Number after Passaging Inactivated Samples






Buffer AVL Results (IFA)

Biodefense solutions to protect our nation

TC 83
Infected
Cells stain
green




Multiplicity of Infection 0.1 IFA Control



Buffer AVL Results (IFA)

Biodefense solutions to protect our nation

**Blue cells
indicate
absence of
infection**






Mock Infection Control



Buffer AVL Results (IFA)

Biodefense solutions to protect our nation






**Green cells indicate
presence of live TC83**

AVL Buffer alone is insufficient to inactivate virus



Buffer AVL Results (IFA)

Biodefense solutions to protect our nation




AVL Buffer with heat sufficiently inactivates virus



SDS Buffer Methodology (JUNV)

Biodefense solutions to protect our nation






SDS Buffer Results (JUNV)

Biodefense solutions to protect our nation

Cell Number for Validation of SDS Inactivation of Junin Infected Cells






SDS Buffer Results (JUNV)

Biodefense solutions to protect our nation


Percent Antibody Positive Cells after Passaging SDS Inactivated Junin Infected Cell Samples





SDS Buffer Results IFA (JUNV)

Biodefense solutions to protect our nation




50 pfu/well JUNV Infection Control




SDS Buffer Results (JUNV)

Biodefense solutions to protect our nation



SDS with 95°C Heat



SDS without 95°C Heat


SDS buffer is effective at inactivating JUNV with or without heat



SDS Buffer Methodology (MERS)

Biodefense solutions to protect our nation


6 well plates+ MOI 5 MERS





SDS Buffer Results (MERS)


Biodefense solutions to protect our nation





SDS Buffer Results (MERS)


Biodefense solutions to protect our nation





SDS Buffer Results IFA (MERS)

Biodefense solutions to protect our nation




MOI 1.5 MERS Infection Control




SDS Buffer Results IFA (MERS)

Biodefense solutions to protect our nation



SDS with 95°C Heat



SDS without 95°C Heat


SDS buffer is effective at inactivating MERS with or without heat




Formalin Methodology (Cells)

Biodefense solutions to protect our nation

6 well plates+ MOI 5 (MERS/EBOV)




Apply 500uL to 145cm² dishes of Veros



6 well plates+ MOI of 1.5 or 0.000005 (MERS/EBOV) in three replicates

Apply 125uL to 145cm² dishes of Veros







Formalin Results (Cells)

Biodefense solutions to protect our nation

Percent Antibody Positive Cells from Formalin Inactivated MERS Infected Cell Monolayer Supernatant Infection




Percent Antibody Positive Cells from Formalin Inactivated Ebola Infected Cell Monolayer Supernatant Infection






Formalin Results (Cells)

Biodefense solutions to protect our nation



Infection from sample of 5 pfu/well MERS at beginning of passaging




Infection from sample of 5 pfu/well EBOV at beginning of passaging




Formalin Results (Cells)

Biodefense solutions to protect our nation



MOI 5 Infected MERS cells following 24 h exposure to formalin





MOI 5 Infected EBOV cells following 24 h exposure to formalin




Formalin Methodology (Tissue)

Biodefense solutions to protect our nation

Liver removed from NHPs succumbing to EBOV infection (titers between 10^7 and 10^8 pfu/mL)




Liver removed from NHPs succumbing to EBOV infection (titers between 10^7 and 10^8 pfu/mL)





Formalin Results (Tissue)

Biodefense solutions to protect our nation






Infection from sample of 5 pfu/well EBOV at
beginning of passaging



Formalin Results (Tissue)

Biodefense solutions to protect our nation




Unfixed, unfiltered EBOV liver samples demonstrating infection



Formalin Results (Tissue)

Biodefense solutions to protect our nation



Fixed, unfiltered EBOV liver samples demonstrating successful inactivation



Trizol LS Methodology (Version 1)

Biodefense solutions to protect our nation

75 μ l NHP Plasma + 25uL of concentrated TC83 in three replicates



Add 300uL Trizol LS



Dilute to 50mL in PBS



Spin/concentrate to ~200-400uL



Spin/wash/concentrate 3x with 12mL PBS each wash



Apply 100% of final contents (~200-400uL) to Vero cells for viral amplification



Incubate 24 hours, passage >50% supernatant, incubated 24 hours, then detect virus by IFA





Trizol LS Methodology (Version 2)

Biodefense solutions to protect our nation

Add TRizol LS treated Zika (R4577T) at 1:1000 in fresh media to a HYPRflask (0.56 ml per flask, 0.14 ml of which is virus and 0.42 ml is TRizol LS)

↓
Incubate 48 hours


↓
Pour off all media and replace with fresh media

↓
Incubate 5 days

↓
Harvest passage one and use 50% for passage two.

↓
Incubate 7 days

↓
Plaque assay analysis





Trizol LS Methodology (Version 3)

Biodefense solutions to protect our nation

Take 4 mL of TRizol treated Zika (R4577T) and use a Zebra column to buffer exchange out the TRizol LS and into completed media (MEM + 10% FBS).

↓
Add the product (100%) of the Buffer exchange into one HYPRflask, ~4mL.

↓
Incubate 7 days

↓
Take 50% of passage 1 for passage 2

↓
Incubate 7 days

↓
Plaque assay analysis






Trizol LS Results (Version 1)

Biodefense solutions to protect our nation

Percent Positive Cells after Passaging Inactivated Samples






Trizol LS Results (Version 1)

Biodefense solutions to protect our nation


Cell Number after Passaging Inactivated Samples





Trizol LS Results (Version 1)

Biodefense solutions to protect our nation




MOI 1.0 TC83 Infection Control



Trizol LS Results (Version 1)

Biodefense solutions to protect our nation



Efficacy of Trizol LS inactivation confirmed across three replicates



Trizol LS Results (Versions 2 and 3)

Biodefense solutions to protect our nation

Sample + 560 mL media	Media Exchanged on Day 2?	Passage 1, Day 7 Observed CPE	Titer on D7 by Plaque Assay	Passage 2, Day 7 Observed CPE	Titer on D7 by Plaque Assay
MEM Only	No	None	None Detected	None	None Detected
Target 100 PFU Zika Actual = 60	No	1+	1.08E+7	4+	1.60E+7
	Yes	1+	5.75E+7	4+	1.65E+7
Target 10 PFU Zika Actual = 4.4	No	Minimal (+/-)	8.00E+7	4+	9.50E+6
	Yes	+/-	3.20E+7	4+	1.13E+7
Target 1 PFU Zika Actual = 0.5	No	None	1.50E+6	4+	5.75E+5
	Yes	None	2.53E+7	4+	3.53E+6
Target 0.1 PFU Zika Actual = 0.05	No	None	None Detected	None	None Detected

Version 2 using dilution method

Version 3 using desalting columns

Sample + 560 mL media	Media Exchanged on Day 2?	Passage 1, Day 7 Observed CPE	Titer on D7 by Plaque Assay (PFU/mL)	Passage 2, Day 7 Observed CPE	Titer on D7 by Plaque Assay (PFU/mL)
Method 1 0.56mL R4577T (Flask A)	Yes	None	None Detected	None	None Detected
Method 1 0.56mL R4577T (Flask B)	Yes	None	None Detected	None	None Detected
Method 1 0.56mL R4577T (Flask C)	Yes	None	None Detected	None	None Detected
Method 2 4mL of R4577T after Zebra buffer Exchange into MEM.	No	None	None Detected	None	None Detected



Conclusions

Biodefense solutions to protect our nation

- USAMRIID has been able to verify efficacy of different toxic chemicals against viruses through dilution, concentration and desalting columns
- Dilution methods are effective for validation of both Trizol LS and formalin based inactivation assays but sample volume is limiting factor
- Both IFA and plaque assays are effective at verifying absence of virus with low LODs (≤ 5 pfu/mL)
- Extensive resources are required at the Institute level to carry out these studies



Acknowledgements

Biodefense solutions to protect our nation

- Cary Retterer
- Tara Kenny
- Rouzbeh Zamani
- Dr. Lou Altamura
- Cindy Rossi
- Brian Kerney
- Jim Jaissle
- Susan Coyne




US Army Medical Research Institute of Infectious Diseases

Mission

Provide leading edge medical capabilities to deter and defend against current and emerging biological threat agents

Vision

To be the leader in the advancement of medical biological defense with world renowned experts dedicated to protecting our military forces and the nation



Medical Biological Defense Insurance Policy for the Nation